**meditec**  
GERMANY . . . . .  
DENTAL / CAD/CAM

# VELVET GRIP

Disposabile NITRILE GLOVES

[www.meditec-germany.de](http://www.meditec-germany.de)  
[www.cerec4smile.de](http://www.cerec4smile.de)  
[www.cad-cam-dental.de](http://www.cad-cam-dental.de)

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## INTRODUCTION CURADEN:

**Curaden AG** ([www.CURADEN.com](http://www.CURADEN.com)) - better health for you - is a public limited company from Switzerland and has been in **existence** for **more than 60 years** (foundation). Our main field of activity is **medical devices (dentists)** and as a result we **would be glad** to cooperate with your company to deliver you premium Nitrile gloves!

## GENERAL INFORMATION:

- **Model:** Nitrile **VELVET GRIP blue** (Premium)
- **Country of Origin:** China
- **Material:** 99,99% pure nitrile
- **Color:** blue (more colors available!)
- **Packaging: GLOVES PACKED FLAT**
  - **Pieces per box:** 100
  - **1 Carton:** 10 boxes (1.000 pcs)
  - **1 pallet:** 120x100cm
    - 1.400 boxes
  - **40 HQ Container**
    - 33.600.500 boxes
- **Sizes:** S, M, L, XL
- **Thickness (material):** min 0.05 mm
- **AQL:** 1.5
- **Dimension (Length x Width):**
  - **Small:** 240 mm x 85 +/- 5mm
  - **Medium:** 240mm x 95 +/- 5mm
  - **Large:** 240mm x 105 +/- 5mm
  - **X-Large:** 240mm x 114 +/- 5mm
- **Features:**
  - **Powder-free:** Yes
  - **Latex-free:** Yes
  - Ambidextrous, textured fingertips
  - Medical grade
  - Non-sterile
- **Storage conditions:** 25C
- **Shelf-life:** 3 years from date of manufacture with the above storage conditions

## CERTIFICATION: (See page 5)

### Manufacturing Accreditation

ISO 9001:2015, ISO 13485:2016

### International Standards

EN 455:1-4; EN 374:1,2,4,5; EN 374-3 (EN 16523-1:2015); EG-Type certification, EN ISO 21420, EN 1186

### Regulation Compliance

CE Marked, FDA 510K

### Declaration of Conformity

EC-REP

## MANUFACTURING ACCREDITATIONS

ISO 9001: 2015

## MANUFACTURING ACCREDITATIONS

### ISO 13485:2016

**Certificate**  
Quality Management System  
EN ISO 13485:2016

Registration No.:  
Organization:

Scope: Design and Development, Manufacture and Distribution of Patient Examination Gloves

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the aforementioned standard are fulfilled. The quality management system is subject to yearly surveillance.

Report No.:  
Effective date: 2021-04-15  
Expiry date: 2024-04-14  
Issue date: 2021-04-13

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**Certificate**  
Quality Management System  
EN ISO 13485:2016

Registration No.:  
Organization:

The scope of certification also covers the following:

No.	Facility	Scope
/01		Distribution of Patient Examination Gloves
/02		Design and Development, Manufacture of Patient Examination Gloves

Report No.:  
Effective date: 2021-04-15  
Expiry date: 2024-04-14  
Issue date: 2021-04-13

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**Certificate**  
Quality Management System  
EN ISO 13485:2016

Registration No.:  
Organization:

The scope of certification also covers the following:

/03	Design and Development, Manufacture of Patient Examination Gloves
/04	Design and Development, Manufacture of Patient Examination Gloves

Report No.:  
Effective date: 2021-04-15  
Expiry date: 2024-04-14  
Issue date: 2021-04-13

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**APPROVED BY GGL**

**Business Direct Products**  
Certification Department

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Application for:  
Certificate No.:  
Registration: EN ISO 13485:2016

Dear Madam or Sir,  
If correct please find the new certificate [attaching file enclosed](#)  
With effective date of the new certificate, the previous certificate becomes invalid.

Best regards,  
Jing Zhang  
TÜV Rheinland  
Certification body

TÜV Rheinland  
LGA Products GmbH  
Tilleystraße 2  
50433 Köln/Lehr  
Germany  
Phone: +49 212 893 2200  
Fax: +49 212 893 2204  
E-Mail: [cert@tuev.com](mailto:cert@tuev.com)  
[www.tuev.com](http://www.tuev.com)  
TÜV Rheinland  
Certification body

**APPROVED BY GGL**

## INTERNATIONAL STANDARDS

### EN 455: 1-4

**Test Report No.**  
dated 20 Oct 2020

ASD Singapore  
For sales  
Inquiry form

**REMARKS:**  
This report is issued solely for the Testing and Certification Registrations of the TUV SUD Group and the General Terms and Conditions of Services of TUV SUD PMS are applicable. In addition, the report is governed by the rules in force within the product.

**SUBJECT:**  
Testing of Gloves submitted by  
on 27 Jul 2020

**TESTED FOR:**

**TEST DATE:**  
22 Jul 2020 to 14 Aug 2020, 20 Oct 2020

**DESCRIPTION OF SAMPLES:**

S/N	Product Description	Colour	Lot No.	Size	Sample received (pieces)	Manufacturer
1	Velvet Disposable Exam Gloves	Blue		S	400	

Lot size as specified by client: 150,000 Pcs 500,000 pieces

**METHOD OF TEST:**

- EN 455-1:2010 Medical gloves for single use Part 1: Requirements and testing for freedom from holes
- EN 455-2:2010 Medical gloves for single use Part 2: Requirements and testing for physical properties
- EN 455-3:2010 Medical glove for single use Part 3: Requirements and testing for biological evaluation Clause 4.4 Freedom from gloves - Clause 4.6 Labeling

Logo: TUV SUD PMS Ltd. 41, Selegie Road, Singapore 347617  
Phone: +65 6336 1100 Fax: +65 6336 1101 Email: info@tuv-sud.com.sg  
Registered Head Office: TUV SUD Asia Pacific Pte Ltd. 100, Cross Street, #02-01, Singapore 110071

**Test Report No.**  
dated 20 Oct 2020

ASD Singapore

**RESULTS:**

Sample: Nitrile Disposable Exam Gloves, L/F No. Size S

**Table 1: Results for EN 455-1:2010**

Clause	Tests	Requirements	Nb. of non-conformities allowed (pieces)	Number tested (pieces)	Actual no. of non-conformities found (pieces)	Inferred results
4	Freedom from holes	Shall not leak	0	310	0	Passed

**Table 2: Results for EN 455-2:2010 Clause 4.5**

Clause	Tests	Requirements (Metric)	Number tested (pieces)	Results (Median)	Inferred results
4	a) Length (mm)	≥ 200	13	246	Passed
	b) Width (mm)	For Size S: 90 ± 5	11	83	Passed
5	a) Strength at break (N)	For nitrile examination glove: ≥ 6.0	13	7.6	Passed
	b) Force at break after challenge testing (N) 7 days at 70±2°C	For nitrile examination glove: ≥ 6.0	11	7.4	Passed

**Table 3: Results for EN 455-3:2010 Clause 7**

Clause	Tests	Requirements	Results	Inferred results
7	Labeling	Manufacturers shall enter the glove expiry date on the packaging with the date of manufacture in accordance with EN ISO 15223-1:2012 and EN 1481:2009+A1:2012. Date of manufacture is defined as the packaging date.	Observed	Passed

**Test Report No.**  
dated 20 Oct 2020

ASD Singapore

**RESULTS (cont'd):**

Sample: Nitrile Disposable Exam Gloves, L/F No. Size S

**Table 4: Results for EN 455-3:2010 Clause 4.4**

Clause	Tests	Requirements	Results / Remarks	Inferred results
4.4	Freedom from gloves	For powder-free gloves: The total quantity of powder residues shall not exceed 2 mg per glove.	0.52 mg per glove	Passed

**Table 5: Results for EN 455-3:2010 Clause 4.6**

Clause	Tests	Requirements	Results
4.6	Labeling	In addition to the labeling specified in EN 12471:2009+A1:2013 and the optional symbols given in EN ISO 15223-1:2012, the following requirements apply.	
		a) manual gloves containing talc shall not be included on the packaging of at least the smallest packaging unit with the EN ISO 15223-1:2012 symbol for talc.	NA
		The labeling shall include the following or equivalent warning statement together with the symbol: (Product) contains natural rubber latex which may cause allergic reactions, including anaphylactic reactions.	NA
		b) the labeling shall include a prominent indication of whether the glove is powder-free or powder.	Comply
		c) sterile powdered gloves shall be labeled with the following or equivalent: CAUTION: Sterile powder shall be removed immediately prior to donning to minimize the risk of airborne latex reactions.	NA
		d) for any medical glove packaging, if used rubber latex (the product labeling shall not include: - any term suggesting relative safety, such as low allergenicity, hypoallergenicity or low pollen; - any unspecified indication of the presence of allergens; - if the manufacturer labels the gloves with the green correct, the process level, measured as specified in 5.3 shall be given.	NA

**Test Report No.**  
dated 20 Oct 2020

ASD Singapore

**REMARKS:**

- Labeling requirements are assessed based on submitted packaging artwork by client on 20 Oct 2020.
- NA: Not applicable for the submitted sample.

**APPENDIX:**

Yeo Peh Tsing Associate Engineer

Lot 08, Yi Engineer Medical Health Services (NAM)

**Photo 1: Nitrile Disposable Exam Gloves, Lot No. Size S**

**Photo 2: Packaging Artwork for Nitrile Disposable**

## INTERNATIONAL STANDARDS

### EN 455: 1-4

**Test Report No.**  
dated 28 Oct 2020

Please note that this Report is issued under the following terms:

- The report applies to the sample of the specific production/shipment given at the time of its homologation. The results are not valid to indicate or imply that they are applicable to other similar units. In addition, such results may not be used to indicate or imply that TÜV SUD PSE approves, recommends or certifies the manufacturer, supplier or user of such production/shipment or that TÜV SUD PSE in any way guarantees the safe performance of the production/shipment unless otherwise stated in this report. No tests were conducted to determine long-term effects of any (bio)compatibility.
- The samples mentioned in this report were submitted by/produced/checked by the Client. TÜV SUD PSE neither assumes, nor is responsible for the accuracy of information on the brand name, model number, origin of manufacture, management or any information reported.
- Nothing in this report shall be interpreted to mean that TÜV SUD PSE has verified or documented any environment or risks from any other testing activity or bodies that may be found in that sample.
- The report shall not be reproduced wholly or in parts and no reference shall be made by the Client to TÜV SUD PSE or to the report or results issued by TÜV SUD PSE in any advertisement or sales promotion.
- Unless otherwise stated, the tests were carried out in TÜV SUD PSE (in Gg) Test Center (in Gg) (Singapore 1002).
- The tests carried out by TÜV SUD PSE and this report are subject to IVD/ISO 13485 Client Terms and Conditions of Business and the Testing and Certification Regulations of the TÜV SUD Group.

Expires 01 September 2023

**Test Report No.**  
dated 28 Oct 2020

**SUBJECT:**  
Testing of Gloves submitted by an EU user 2020.

**TESTED FOR:**

**TEST DATE:**  
22 Jul 2020 to 14 Aug 2020, 28 Oct 2020

**DESCRIPTION OF SAMPLES:**

S/N	Product Description	Colour	Lot No.	Size	Sample quantity (pieces)	Manufacturer
1	White Disposable Exam Gloves	White	M	M	400	

Lot size as specified by client: 100,000 to 100,000 gloves

**METHOD OF TEST:**

- EN 455-1:2010 Medical gloves for single use Part 1: Requirements and testing for structures from tables
- EN 455-2:2010 Medical gloves for single use Part 2: Requirements and testing for physical properties
- EN 455-3:2010 Medical glove for single use Part 3: Requirements and testing for biological evaluation
  - Clause 4.4 Powder-free gloves
  - Clause 4.6 Labelling

**Test Report No.**  
dated 28 Oct 2020

**RESULTS**

Sample: White Disposable Exam Gloves.

**Table 1: Results for EN 455-1:2010**

Clause	Test	Requirements	No. of Parameters tested (pieces)	Number failed (pieces)	Actual no. of non-compliant (failed) pieces	Passed results.
4	Friction from latex	Shall not leak	16	310	3	Passed

**Table 2: Results for EN 455-2:2010 Clauses 4.3**

Clause	Test	Requirements (Method)	Number tested (pieces)	Results (Median)	Intended results
4	Dimensions (a) Length (mm)	≥ 240	13	240	Passed
	(b) Width (mm)	For 2000 M 90 x 10	10	93	Passed
	Strength at Force at break (N)	For nitrile non-sterile gloves: ≥ 5.0	15	8.8	Passed
5	(c) Force at break after challenge testing (N) 7 days at 23±0.5°C	For nitrile examination gloves: ≥ 5.0	13	7.1	Passed

**Table 3: Results for EN 455-3:2010 Clause 3**

Clause	Test	Requirements	Results	Intended results.
1	Labelling	Manufacturers shall label the glove and/or the packaging with the date of manufacture in accordance with EN ISO 15223-1:2013 and EN ISO 15424:2004+A1:2013. Date of manufacture is defined as the packaging date.	Observed	Passed

**Test Report No.**  
dated 28 Oct 2020

**RESULTS**

Sample: White Disposable Exam Gloves, Lot No.

**Table 4: Results for EN 455-3:2010 Clause 4.4**

Clause	Test	Requirements	Results / Remarks	Intended results
4.4	Powder-free gloves	For powder-free gloves, the total quantity of residual powders shall not exceed 2 mg per glove.	1.08 mg per glove	Passed

**Table 5: Results for EN 455-3:2010 Clause 4.6**

Clause	Test	Requirements	Results
4.6	Labelling	In addition to the labelling specified in EN ISO 15424:2004+A1:2013 and the relevant symbols given in EN ISO 15223-1:2013, the following requirements apply: a) medical gloves containing natural rubber latex shall be labelled on the packaging of at least the smallest packaging unit with the EN ISO 15223-1:2013 symbol for latex. The labelling shall include the following or equivalent warning statement together with the symbol: "Products containing natural rubber latex which may cause allergic reactions; industry employment suggested." b) the labelling shall include a prominent indication of whether the glove is powdered or powder-free. c) nitrile examination gloves shall be labelled with the following or equivalent: "CAUTION: Surface powder shall be removed" specifically prior to undertaking operative procedures in order to maintain the risk of allergic disease reactions. d) for any medical glove containing natural rubber latex for product labelling shall not include: - any term suggesting relative safety, such as the statements: "hypoallergenic" or "low protein"; - any unqualified indication of the presence of latex; or e) if the manufacturer labels the gloves with the generic content, the process limit, measured as specified in 2.3.6.40, is given. Intended results: Passed	Comply

## INTERNATIONAL STANDARDS

### EN 455: 1-4

**Test Report No.**  
dated 28 Oct 2020

Please note that this Report is issued under the following terms:

- The report applies to the sample of the specific production/shipment given at the time of its homologation. The results are not valid to indicate or imply that they are applicable to other similar units. In addition, such results may not be used to indicate or imply that TÜV SUD PSE approves, recommends or certifies the manufacturer, supplier or user of such production/shipment or that TÜV SUD PSE in any way guarantees the safe performance of the production/shipment unless otherwise stated in this report. No tests were conducted to determine long-term effects of any (bio)compatibility.
- The samples mentioned in this report were submitted by/produced/checked by the Client. TÜV SUD PSE neither assumes, nor is responsible for the accuracy of information on the brand name, model number, origin of manufacture, management or any information reported.
- Nothing in this report shall be interpreted to mean that TÜV SUD PSE has verified or documented any environment or risks from any other testing activity or bodies that may be found in that sample.
- The report shall not be reproduced wholly or in parts and no reference shall be made by the Client to TÜV SUD PSE or to the report in results forwarded by TÜV SUD PSE in any advertisement or sales promotion.
- Unless otherwise stated, the tests were carried out in TÜV SUD PSE (in Gg) Test Center (in Gg) (Singapore 1002).
- The tests carried out by TÜV SUD PSE and this report are subject to IVD/ISO 13485 Client Terms and Conditions of Business and the Testing and Certification Regulations of the TÜV SUD Group.

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**SUBJECT:**  
Testing of Gloves submitted by an EU user 2020.

**TESTED FOR:**

**TEST DATE:**  
22 Jul 2020 to 14 Aug 2020, 28 Oct 2020

**DESCRIPTION OF SAMPLES:**

S/N	Product Description	Colour	Lot No.	Size	Sample quantity (pieces)	Manufacturer
1	White Disposable Exam Gloves	White	M	M	400	

Lot size as specified by client: 100,000 to 100,000 gloves

**METHOD OF TEST:**

- EN 455-1:2010 Medical gloves for single use Part 1: Requirements and testing for structures from latex
- EN 455-2:2010 Medical gloves for single use Part 2: Requirements and testing for physical properties
- EN 455-3:2010 Medical glove for single use Part 3: Requirements and testing for biological evaluation
  - Clause 4.4 Powder-free gloves
  - Clause 4.6 Labelling

**Test Report No.**  
dated 28 Oct 2020

**RESULTS**

Sample: White Disposable Exam Gloves

**Table 1: Results for EN 455-1:2010**

Clause	Test	Requirements	No. of Parameters tested (pieces)	Number failed (pieces)	Actual no. of non-compliant (failed) pieces	Passed results
4	Friction from latex	Shall not fail	16	316	3	Passed

**Table 2: Results for EN 455-2:2010 Clauses 4.3**

Clause	Test	Requirements (Method)	Number tested (pieces)	Results (Median)	Intended results
4	Dimensions (a) Length (mm)	≥ 240	15	240	Passed
	(b) Width (mm)	For 2000 M 88 x 10	15	93	Passed
	Strength at Force at break (N)	For nitrile non-sterile gloves: ≥ 5.0	15	8.8	Passed
5	(c) Force at break after challenge testing (N) 7 days at 23±0.5°C	For nitrile examination gloves: ≥ 5.0	15	7.1	Passed

**Table 3: Results for EN 455-3:2010 Clause 3**

Clause	Test	Requirements	Results	Intended results
1	Labelling	Manufacturers shall label the glove and/or the packaging with the date of manufacture in accordance with EN ISO 15223-1:2013 and EN ISO 15424:2004+A1:2013. Date of manufacture is defined as the packaging date.	Observed	Passed

**Test Report No.**  
dated 28 Oct 2020

**RESULTS**


Sample: White Disposable Exam Gloves, Lot No.

**Table 4: Results for EN 455-3:2010 Clause 4.4**

Clause	Test	Requirements	Results / Remarks	Intended results
4.4	Powder-free gloves	For powder-free gloves, the total quantity of residual powders shall not exceed 2 mg per glove.	1.08 mg per glove	Passed

**Table 5: Results for EN 455-3:2010 Clause 4.6**

Clause	Test	Requirements	Results
4.6	Labelling	In addition to the labelling specified in EN ISO 15424:2004+A1:2013 and the relevant symbols given in EN ISO 15223-1:2013, the following requirements apply: a) medical gloves containing natural rubber latex shall be labelled on the packaging of at least the smallest packaging unit with the EN ISO 15223-1:2013 symbol for latex. The labelling shall include the following or equivalent warning statement together with the symbol: "Product contains natural rubber latex which may cause allergic reactions; industry employment suggested." b) the labelling shall include a prominent indication of whether the glove is powdered or powder-free. c) nitrile examination gloves shall be labelled with the following or equivalent: "CAUTION: Surface powder; shall be removed" specifically prior to undertaking operative procedures in order to maintain the risk of allergic disease reactions. d) for any medical glove containing natural rubber latex for product labelling shall not include: - any term suggesting relative safety, such as the statements: "hypoallergenic" or "low protein"; - any unqualified indication of the presence of latex; or e) if the manufacturer labels the gloves with the generic content, the process limit, measured as specified in 5.2.3.040, is given. Intended results: Passed	NA NA Comply NA NA Passed

## INTERNATIONAL STANDARDS

### EN 455:1-4

**Test Report No.**  
dated 20 Oct 2020

**RESULTS (cont'd)**  
Sample: Nitrile Disposable Exam Gloves, Lot No.

**Table 6: Results for EN 455-2:2013 Clause 4.4**

Clause	Tests	Requirements	Results / Remarks	Interim Results
4.4	Powder-free gloves	For powder-free gloves, the total quantity of powder residues shall not exceed 0.1 mg per glove.	0.08 mg per glove	Passed

**Table 8: Results for EN 455-3:2013 Clause 4.6**

Clause	Tests	Requirements	Results
4.6	Labeling	In addition to the labeling specified in EN ISO 13998-4:2013 and the relevant standards given in EN ISO 15223-1:2012, the following requirements apply: a) medical gloves containing natural rubber latex shall be labeled on the packaging of at least the smallest packaging unit with the EN ISO 15223-1:2012 symbol for latex. The labeling shall include the following or equivalent warning statement together with the symbol: (Product) contains natural rubber latex which may cause allergic reactions including anaphylaxis to latex. b) The labeling shall include a statement indicating whether the glove is powder-free or powder-free. c) sterile powdered gloves shall be labeled with the following or equivalent: "CAUTION: Surface powder shall be removed sterily prior to undertaking operative procedures in order to minimize the risk of allergic disease reactions." d) for any reusable powder-free rubber gloves the product labeling shall not include: - any text regarding re-use safety (such as one-way, hydrophobicity or the gloves), - any unqualified indication of the presence of chlorine, e) if the manufacturer labels the gloves with the product content, the product unit, measured as specified in 5.1 shall be given.	Passed

**Test Report No.**  
dated 20 Oct 2020

**REMARKS:**

- Labeling requirements are assessed based on submitted packaging artwork by client on 20 Oct 2020.
- NA: Not applicable for the submitted sample.

Pro Peh Kang Associate Engineer  
 Eng Chyi Engineer  
 Medical Health Services (MMS)

**APPENDIX:**

Photo 1: Nitrile Disposable Exam Gloves, Lot No. Blue

Photo 2: Packaging Artwork for Nitrile Disposable Exam Gloves, Lot No.

**Test Report No.**  
dated 20 Oct 2020

**Please note that this Report is issued under the following terms:**

- This report applies to the sample of the specific product/glove given at the time of its submission. The results are not valid to other samples that may be submitted at other times. In addition, such results may not be used to release or stop the TUV SUD PSE approval, requirements or evidence the manufacturer, supplier or user of such product/glove, or stop TUV SUD PSE in any way, regardless the size performance of the product/glove. Unless otherwise stated in this report, no tests were conducted to determine long-term effects of using the specific product/glove.
- The complete responsibility of this report shall remain with the manufacturer or the client. TUV SUD PSE assumes no responsibility for the accuracy of information on its label name, model number, size of production, origin or any other information supplied.
- Nothing in this report shall be interpreted to mean that TUV SUD PSE has verified or ascertained any requirement or made any other finding as a party or before that may be found on the sample.
- This report shall not be reproduced, wholly or in parts, and its address shall be made by the Client to TUV SUD PSE in the report or results furnished by TUV SUD PSE in any subsequent or similar production.
- Unless otherwise stated, this has been carried out by TUV SUD PSE Pte Ltd, No. 1 Science Park Drive, Singapore 119221.
- The tests carried out by TUV SUD PSE and this report are subject to TUV SUD PSE's General Terms and Conditions of Business and the Testing and Certification Regulations of the TUV SUD Group.

Effective 11 September 2020

**Test Report No.**  
dated 20 Oct 2020

**NOTE:** This report is issued subject to the Testing and Certification Regulations of the TUV SUD Group and the General Terms and Conditions of Business of TUV SUD PSE Pte Ltd. In addition, this report is governed by the Intellectual Property Rights.

**SUBJECT:**  
Testing of Gloves submitted by client on 27 Jul 2020.

**TESTED FOR:**

**TEST DATE:**  
22 Jul 2020 to 14 Aug 2020; 20 Oct 2020

**DESCRIPTION OF SAMPLES:**

SN	Product Description	Color	Lot No.	Size	Sample received (Quantity)	Manufacturer
1	Nitrile Disposable Exam Gloves	Blue	NA	401	401	

Lot size as specified by client: 130,000 to 500,000 pieces.

**METHOD OF TEST:**

- EN 455-1:2010 Assestle gloves for single use  
Part 1: Requirements and testing for treatment from latex
- EN 455-2:2013 Medical gloves for single use  
Part 2: Requirements and testing for physical properties
- EN 455-3:2013 Medical glove for single use  
Part 3: Requirements and testing for biological evaluation
  - Clause 4.4 Powder-free gloves
  - Clause 4.6 Labeling

## INTERNATIONAL STANDARDS

### EN 455:1-4

**Test Report No.**  
dated 20 Oct 2020

**RESULTS**

Sample: Nitrile Disposable Exam Glove, Lot No. X1

**Table 1: Results for EN 455:1:2012**

Clause	Tests	Requirements	No. of non-complies allowed (pieces)	Number tested (pieces)	Actual no. of non-complies found (pieces)	Inferred results
4	Freedom from holes	Shall not leak	10	315	2	Passed

**Table 2: Results for EN 455:2:2013 Clauses 4-6**

Clause	Tests	Requirements (Method)	Number tested (pieces)	Results (Meets)	Inferred results
4	a) Length (mm)	4.340	55	248	Passed
		b) Width (mm)	For size M: 8.5 to 9.5	114	
5	a) Tensile strength at break (N)	For nitrile examination gloves: 8.0	53	8.8	Passed
		b) Tensile strength at break after challenge (N) 7 days at 120% TC	8.0	7.0	

**Table 3: Results for EN 455:3:2013 Clause 7**

Clause	Tests	Requirements	Results	Inferred results
7	Labelling	Manufacturers shall ensure the glove inside the packaging with the date of manufacture in accordance with EN ISO 15224-1:2017 and EN ISO 20014:2013. Date of manufacture is defined as the packaging date.	Observed	Passed

**Test Report No.**  
dated 20 Oct 2020

**REMARKS (cont'd):**

Sample: Nitrile Disposable Exam Glove, Lot No. X1

**Table 4: Results for EN 455:3:2013 Clause 8.1**

Clause	Tests	Requirements	Results (Remarks)	Inferred results
8.1	Protein-free gloves	For powder-free gloves: The total quantity of protein residues shall not exceed 2 mg per glove	0.61 mg per glove	Passed

**Table 5: Results for EN 455:3:2013 Clause 8.2**

Clause	Tests	Requirements	Results	
8.2	Labelling	In addition to the labelling specified in EN ISO 20014:2013 and the relevant symbols given in EN ISO 15225-1:2012, the following requirements apply:	Inferred results	
		a) The nitrile glove containing natural rubber latex shall be indicated on the packaging or at least the sealed packaging with the EN ISO 15225-1:2012 symbol for latex.		NA
		b) The labelling shall include the following or equivalent warning information together with the symbol: (Products) contain latex rubber latex which may cause allergic reactions, including anaphylactic responses		NA
		c) The labelling shall include a prominent indication of whether the glove is powdered or powder-free		Comply
		d) Specific powdered glove label be labelled with the following in maximum CAUTION: Surface powder may be removed especially prior to undertaking operative procedures in order to minimize the risk of adverse tissue reactions		NA
		e) For any nitrile glove containing natural rubber latex the product labelling shall not include: <ul style="list-style-type: none"> <li>- any false suggestive relative safety, such as low temperature, hypoallergenicity or skin protectant</li> <li>- any ambiguous indication of the presence of allergen</li> </ul>		NA
f) If the manufacturer labels the gloves with the protein content, the protein limit, measured as specified in 8.1 shall be given.	NA			

**Test Report No.**  
dated 20 Oct 2020

**REMARKS:**

- Labelling requirements are assessed based on submitted packaging artwork by client on 20 Oct 2020.
- NA: Not applicable for the submitted sample.

**APPENDIX:**

Photo 1: Nitrile Disposable Exam Glove, Lot No. X1

Photo 2: Packaging Artwork for Nitrile Disposable Exam Glove, Lot No. X1

**Test Report No.**  
dated 20 Oct 2020

Please note that this Report is issued under the following terms:

- This report applies to the specific health professional's product as the object of its investigation. This report will not be used to indicate or imply that they are suitable for other uses. In addition, such results may not be used to indicate or imply that TUV SUD PSE approves, recommends or endorses the manufacturer, supplier or user of said product(s) or that TUV SUD PSE is any way guaranteeing the safe performance of the product(s). Unless otherwise stated in this report, no item was produced to determine long-term effects of using the specific product(s).
- The samples considered in this report were submitted to the Client by TUV SUD PSE. Further samples are recommended, for the accuracy of information on the stated items, include samples, origin of production, composition or any alternative material.
- Nothing in this report shall be construed to mean that TUV SUD PSE has verified or determined any statement or made from specific testing activities to follow the use or test of the sample.
- This report will not be reproduced, stored in a retrieval system or otherwise disseminated to the Client by TUV SUD PSE or to the report of health professionals in any other form or by any other means.
- Unless otherwise stated, the tests were carried out in TUV SUD PSE (TU) Ltd, Schwanen Drive, Singapore 110271.
- The work carried out by TUV SUD PSE and the user acknowledged by TUV SUD PSE (TU) Ltd and the Client of the Terms and Conditions of Business and the Testing and Certification Regulations apply to this Report.

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## INTERNATIONAL STANDARDS

### EN 374 1,5

Produkt Product			TÜVRheinland®	
Prüfbericht-Nr.: Test Report No.		Seite 21 von 24 Page 21 of 24		
Abkürz.	EN ISO 374-1:2016 + A1:2018	Messergebnisse - Bemerkungen	Bewertung	
Clause	Anforderungen - Prüfungen / Requirements - Tests	Measuring results - Remarks	Evaluation	
	The following minimum information shall be supplied - Name and full address of manufacturer or his authorized representative			
	- Glove designation			
	- Information on available size range	given		
	- Reference to the relevant specific European standard, paragraph with performance level	given		
	- If the expected shelf life of the gloves is reduced by aging, the expiration date have to be added or information regarding shelf life	given 3 years		
	- If protection is only given for part of gloves, information have to be added	not		
	- possible problems	given		
	- instruction for use for gloves and also for use with combination of other PPE	NA		
	- A list of the substances contained in the glove which are known to cause allergies	NA		
	- care symbols or explanations	not		
	- storage instructions	given		
	- Type of packaging suitable for transport, if relevant	given		
	- Name and address of the testing laboratory and/or its number	given or not		
	Furthermore, a basic explanation shall be given to assist comprehension of the relevant performance levels, and the standard(s) to which they refer shall be included	given		

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Produkt Product			TÜVRheinland®	
Prüfbericht-Nr.: Test Report No.		Seite 23 von 24 Page 23 of 24		
Abkürz.	EN ISO 374-1:2016 + A1:2018	Messergebnisse - Bemerkungen	Bewertung	
Clause	Anforderungen - Prüfungen / Requirements - Tests	Measuring results - Remarks	Evaluation	
	Following information shall be added in user instructions - the results of Penetration, Degradation, Permeation - the list of all the chemicals to which the protective gloves have been tested - the performance levels obtained in permeation testing	given		
	The following warnings shall be added in the user instructions: - This information does not reflect the actual duration of protection in the workplace and the differentiation between industries and pure chemical	given		
	- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only checked in cases where the glove is equal to or over 400 mm - where the cuff is tested later and relate only to the chemical tested. It can be different if the chemical is used in a mixture.	given		
	- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.	given		
	- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.	given		
	- Before usage, inspect the gloves for any defect or imperfections.	given		
	For reusable gloves, the manufacturer shall provide the relevant instructions for decontamination. If there is no information about decontamination, then it is intended for single use only and the following warning shall be added: "For single use only."	NA		

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## INTERNATIONAL STANDARDS

### EN 374 1,5

Produkt Product			TÜVRheinland®		
Prüfbericht-Nr.: Test Report No.		Seite 21 von 24 Page 21 of 24			
Ansatz	EN ISO 374-1:2016 + A1:2018	Messergebnisse - Bemerkungen	Bewertung		
Clause	Anforderungen - Prüfungen / Requirements - Tests	Measuring results - Remarks	Evaluation		
	The following minimum information shall be supplied - Name and full address of manufacturer or his authorized representative				
	- Glove designation				
	- Information on available size range	given			
	- Reference to the relevant specific European standard, paragraph with performance level	given			
	- If the expected shelf life of the gloves is reduced by aging, the expiration date have to be added or information regarding shelf life	given 3 years			
	- If protection is only given for part of gloves, information have to be added	not			
	- possible problems	given			
	- instruction for use for gloves and also for use with combination of other PPE	NA			
	- A list of the substances contained in the glove which are known to cause allergies	NA			
	- care symbols or explanations	not			
	- storage instructions	given			
	- Type of packaging suitable for transport, if relevant	given			
	- Name and address of the testing laboratory and/or its number	given or not			
	Furthermore, a basic explanation shall be given to assist comprehension of the relevant performance levels, and the standard(s) to which they refer shall be included	given			

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Produkt Product			TÜVRheinland®		
Prüfbericht-Nr.: Test Report No.		Seite 23 von 24 Page 23 of 24			
Ansatz	EN ISO 374-1:2016 + A1:2018	Messergebnisse - Bemerkungen	Bewertung		
Clause	Anforderungen - Prüfungen / Requirements - Tests	Measuring results - Remarks	Evaluation		
	Following information shall be added in user instructions - the results of Penetration, Degradation, Permeation - the list of all the chemicals to which the protective gloves have been tested - the performance levels obtained in permeation testing	given			
	The following warnings shall be added in the user instructions: - This information does not reflect the actual duration of protection in the workplace and the differentiation between industries and pure chemical	given			
	- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only checked in cases where the glove is equal to or over 400 mm - where the cuff is tested later and relate only to the chemical tested. It can be different if the chemical is used in a mixture.	given			
	- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.	given			
	- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.	given			
	- Before usage, inspect the gloves for any defect or imperfections.	given			
	For reusable gloves, the manufacturer shall provide the relevant instructions for decontamination. If there is no information about decontamination, then it is intended for single use only and the following warning shall be added: "For single use only."	NA			

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Customer details: SATRA reference: Your reference: Date of report: 23 March 2021 Samples received: 1 March 2021 Date(s) work carried out: 4-23 March 2021

### TECHNICAL REPORT

Subject: EN ISO 21420: 2020 size & clarity & innocuousness test, EN ISO 374-2: 2019 air leak and water leak, EN ISO 374-5: 2016 viruses on Disposable Powder Free Nitrile Gloves referenced as colour: blue, size: X(S)-6, S(5)-7, M(7)-8, L(9)-9, XL(9)-10.

**Conditions of Issue:**  
This report may be forwarded to other parties provided that it is not charged in any way. It must not be published, for example by including it in advertisements, without the prior written permission of SATRA.  
Results given in this report refer only to the samples submitted for analysis and tested by SATRA. Comments are for guidance only.  
A satisfactory test report in no way implies that the product tested is approved by SATRA and no warranty is given as to the performance of the product tested. SATRA shall not be liable for any subsequent loss or damage incurred by the client as a result of information supplied in this report.  
The uncertainty of the results (U<sub>M</sub>) in this report is based on a standard uncertainty multiplied by a coverage factor (k=2) which provides a coverage probability of approximately 95%.

Report signed by: Gladys He  
Position: Technologist  
Department: China Testing

(Page 1 of 5)

**SATRA TECHNOLOGY**

### TECHNICAL REPORT

**WORK REQUESTED**  
Samples described as Disposable Powder Free Nitrile Gloves referenced as colour: blue, size: X(S)-6, S(5)-7, M(7)-8, L(9)-9, XL(9)-10 were received by SATRA on 1 March 2021 for testing in accordance with EN ISO 21420: 2020, EN ISO 374-2: 2019 and EN ISO 374-5: 2016.

**SAMPLE SUBMITTED**

Samples described as Disposable Powder Free Nitrile Gloves referenced as colour: blue

**TESTING REQUESTED**  
EN ISO 21420: 2020 Clause 5.1 – Sizing and measurement of gloves  
EN ISO 21420: 2020 Clause 5.2 – Clarity  
EN ISO 374-2: 2019 Clause 7.2 – Air leak  
EN ISO 374-2: 2019 Clause 7.3 – Water leak  
EN ISO 374-5: 2016 Clause 5.3 – Protection against viruses (ISO 19904: 2004 Procedure B)  
EN ISO 21420: 2020 Clause 4.2 – Innocuousness of protective gloves

**CONCLUSION**  
The samples described as Disposable Powder Free Nitrile Gloves referenced as colour: blue, size: X(S)-6, S(5)-7, M(7)-8, L(9)-9, XL(9)-10 were found to achieve the following results:  
EN ISO 21420: 2020 Clause 5.1 – See below table  
EN ISO 21420: 2020 Clause 5.2 – Level 5  
EN ISO 374-2: 2019 Clause 7.2 – Pass  
EN ISO 374-2: 2019 Clause 7.3 – Pass  
EN ISO 374-5: 2016 Clause 5.3 – Pass  
EN ISO 21420: 2020 Clause 4.2 – Pass PAVs, pH value and DMF  
Detailed results are included on the following page(s)

Date: 23 March 2021 (Page 2 of 5)

Signature: Gladys He  
China Testing

**SATRA TECHNOLOGY**

### TECHNICAL REPORT

**Testing**  
Testing was carried out in accordance with EN ISO 21420:2020 and EN ISO 374-2: 2019  
Samples for testing were conditioned for at least 24 hours in a conditioned environment maintained at (23±2) °C and (50±5) % relative humidity.

**Requirements**  
Table 1 – Requirements for EN ISO 21420: 2020 Clause 5.2 Clarity

Performance level	1	2	3	4	5
Diameter of clarity pin (mm)	11.0	9.5	8.0	6.5	5.0

Table 2 – Requirements for EN ISO 374-2: 2019

Clause 7.2 Air leak	No leak to be detected
Clause 7.3 Water leak	No leak to be detected

**Test Results**  
Table 3 – EN ISO 21420: 2020 Test Results

Clause / Test	Requirement	Test Results	U <sub>M</sub> (See note #)	Result				
5.1 Glove length, comfort and fit	N/A	Size	Length (mm)	± 1.10 mm	N/A			
		1	2			3		
		5-6	237			234	240	
		Comfortable on fit	6-7			230	235	236
		Comfortable on fit	7-8			248	245	240
		Comfortable on fit	8-9			237	235	240
		Comfortable on fit	9-10			254	258	260
5.2 Clarity	See table 1	Size	Minimum pin diameter / mm	N/A				
		5-6	5.0					
		6-7	5.5					
		7-8	5.0					
		8-9	5.0					

Date: 23 March 2021 (Page 3 of 5)

Signature: Gladys He  
China Testing

**SATRA TECHNOLOGY**

### TECHNICAL REPORT

Table 4 – EN ISO 374-2: 2019 Test Results

Clause / Test	Test Results	Unit (See note #)	Result
7.2 Air leak test	Total air pressure used	3.1 kPa	N/A
	Sample size	Leak(s)	
	5-6	No leaks detected	
	6-7	No leaks detected	
	7-8	No leaks detected	
7.3 Water leak test	Sample size	Leak(s)	N/A
	5-6	No leaks detected	
	6-7	No leaks detected	
	7-8	No leaks detected	
	8-9	No leaks detected	

**Additional Information / Notes**  
Note # – Estimated uncertainty of measurement applied at point of test (e.g. to applied force or to tolerance limit) to ensure product meets requirements of the standard

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Determine the degradation for each of the three glove specimens against each specific chemical or chemical mixture.

Determine the standard deviation (SD) of the degradation for the three gloves.

Any changes such as swelling, shrinking, softening, deformation, enlargement, softening, melting, expansion, expansion, hardening/bleaching, stretching shall be noted and described for information.

P2 Chemikalie / chemical  
Reinheit / purity: 40%  
D01 -12,5  
D02 -11,9  
D03 -11,9  
DR -4,3  
SD 14,6

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## MANUFACTURING ACCREDITATIONS

CE 0197

**Business Stream Products**  
Certification Department

**TÜV Rheinland**  
LGA  
Productivity Right

100 Pfeiderer-Str. Industrial Estate - Bielefeld, Germany

Address:  
Tel. +49 51 900 0200  
Fax +49 51 900 0201  
Internet: www.tuev.com

Application for: EU type-examination certificate PPE  
Certificate No.:  
Device: Protective gloves against chemicals and microorganisms according to EN ISO 374-1+A1:2018 and EN ISO 374-5:2018  
Type:  
Test requirement: (EN)Reg 426/2016  
EN ISO 374-1:2018+A1  
EN ISO 374-5:2018

Dear Madam or Sir,  
A specimen of above mentioned product has been tested and found to be technically in compliance with § 20 and § 21 of the ProdSG (in German Product Safety Law).  
The certificate is issued with the reservation that the licensee holder applies all information required in § 6 of the ProdSG related to name of the manufacturer and, if need be, its authorized representative / the importer including their respective address on the product, its packaging and/or the user's manual prior to marketing the product in the European Economic Area.  
Enclosed please find the certificate of approval No. \_\_\_\_\_  
Kind regards  
Certification body  
C. Albrecht  
Dipl.-Ing. C. Albrecht  
Test sample no. documentation available

TÜV-Rheinland  
LGA-Produkt-Service  
Telefon +49 51 900 0200  
Fax +49 51 900 0201  
Web: www.tuev.com  
E-Mail: ggg@tuev.com  
Service of Management  
TUV-RG  
Key-Market: Industries  
TUV-RG  
Dipl.-Ing. Albrecht  
Chairman of the Supervisory Board  
TUV-RG  
Key-Market: Industries  
Dipl.-Ing. Albrecht

**CERTIFICATE**  
EU Type-Examination Certificate  
Regulation 2016/425/EU  
Personal Protective Equipment

**TÜV Rheinland**

Registration No.:  
Report No.:

Holder:

Product: Protective gloves against chemicals and microorganisms according to EN ISO 374-1+A1:2018 and EN ISO 374-5:2018

Identification: disposable gloves  
Colour: blue  
Material: nitrile, wall thickness 0,55 mm  
Size: S (4,5) - XL (9)  
Performance parameter: Type C, class 4, S: 3x02 419  
- PPE Category III - obligatory monitoring module C2 -  
The EU type-examination certificate refers to the above mentioned product. This is to certify that the product complies with the essential requirements of Annex II of the regulation 2016/425/EU. This certificate does not imply assessment of the production of the product and does not prevent the use of a TÜV Rheinland mark of conformity. The holder is entitled to use this certificate in connection with the declaration of conformity in accordance with Annex IX.

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CE The CE marking may be used if all relevant and effective EC Directives are complied with.

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**CLIENT NAME**  
**CLIENT ADDRESS**

**TEST PERIOD** 01-Mar-2018-08-Mar-2019

**TEST REQUEST** In accordance with Council of Europe Res AP (2004) 4

**CONCLUSION** **PASS**  
The submitted sample was found to comply with the overall migration requirements as stated in European Regulation Res AP (2004) 4 on rubber to be used for food contact applications.

Prepared By: *Cynthia Cao*  
(Cynthia Cao)  
Report Draftsman

Authorized By: *Lei*  
(Lei)  
Authorized Signatory

**APPROVED BY**  
GGL

Notes: (1) Shows "Limit & Condition as mentioned" (2) This result is only for the items tested (3) The test method is standardized method in full unless the written approval of the laboratory (4) Through the agreement of the laboratory, the user may choose to use the test results for other applications

Client/In-house Laboratory: TÜV SÜD Products Testing (Shanghai) Co., Ltd. 8-204, No. 1008 Da Hua Road, Minhang District Shanghai 201108 P.R. China  
Phone: +86 (21) 597 675 Fax: +86 (21) 597 676 Email: test@tuv-sud.com.cn Website: www.tuv-sud.com.cn  
Regional Head Office: TÜV SÜD Certification and Testing (Shanghai) Co., Ltd. No. 20, Jinqiang Road (Shanghai) 200070, P.R. China  
TUV

Test Report No.:  
Report Date: 11 April 2019

**RECEIPT DATE / TEST DATE**  
01-Mar-2019 / 01-Mar-2019

**THE FOLLOWING SAMPLES WERE SUBMITTED BY/ ON BEHALF OF THE CLIENTS AS**

Sample Name: Powder free nitrile glove, blue  
Sample Specification: Medium  
Batch/Article: Manufacture

SAMPLE NO.	DESCRIPTION	PHOTOGRAPH
	Blue glove	

**TEST METHOD(S)**

1. For material: Rubber  
- Overall migration test for compliance with European Regulation Res AP (2004) 4 on rubber to be used for food contact applications.  
- As specified in REGULATION (EU) No 1005/11 and its amendments, with reference to EN 1186, Part 2 (Test methods for overall migration into aqueous food simulants) / EN 1186, Part 3 (Test methods for overall migration into aqueous food simulants by total ion chromatography) / EN 1186, Part 14 (substitute test)

**TEST RESULT(S)**

1. Overall Migration Test - with reference to EN 1186: Part 2, Part 3 & Part 14

Simulated Food	Test Condition	Result (mg/kg)	Maximum Permissible Limit (mg/kg)
10% Ethanol	75°C for 2 hours	<= 0.0	60
2% Acetic acid	75°C for 2 hours	<= 0.0	60
20% Ethanol	80°C for 2 hours	<= 0.0	60
Isopropyl	85°C for 0.5 hour	<= 0.0	60

Note: 1. mg/kg denotes milligram per kilogram foodstuff.  
2. Specification is quoted from European Regulation Res AP (2004) 4 on rubber to be used for food contact applications.  
3. \* denotes less than

Note: This report is for internal use only by the client.

**END OF THE TEST REPORT**

Client/In-house Laboratory: TÜV SÜD Products Testing (Shanghai) Co., Ltd. 8-204, No. 1008 Da Hua Road, Minhang District Shanghai 201108 P.R. China  
Phone: +86 (21) 597 675 Fax: +86 (21) 597 676 Email: test@tuv-sud.com.cn Website: www.tuv-sud.com.cn  
Regional Head Office: TÜV SÜD Certification and Testing (Shanghai) Co., Ltd. No. 20, Jinqiang Road (Shanghai) 200070, P.R. China  
TUV

**APPROVED BY**  
GGL

## REGULATION COMPLIANCE

### FDA 510K

The screenshot shows the FDA's 510(k) Premarket Notification database entry for 'Polymer Patient Examination Glove'. The page includes a search bar, navigation tabs, and a detailed table of product and regulatory information.

New Search		Back To Search Results
Device Classification Name	Polymer Patient Examination Glove	
510(k) Number		
Device Name	Powder Free Nitrile Patient Examination Glove, Blue Colored, Non Sterile, Tested For Use With Chemotherapy Drugs	
Applicant		
Applicant Contact Correspondent		
Correspondent Contact		
Regulation Number		
Classification Product Code	L26	
Subsequent Product Code	L26	
Date Received	06/23/2017	
Decision Date	11/15/2017	
Decision	Substantially Equivalent (SESE)	
Regulation Medical Specialty	General Hospital	
510k Reviewer Panel	General Hospital	
Summary	Summary	
Type	Traditional	
Reviewed By Third Party	No	
Combination Product	No	

Page Last Updated: 06/28/2021  
 Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.  
 Language Assistance Available: Español | 繁體中文 | Tiếng Việt | বাংলা | Tagalog | ਪੰਜਾਬੀ | العربية | Kreyòl Ayisyen | Français | Polski | Português | हिन्दी | Deutsch | 日本語 | العربية | English

## DECLARATION OF CONFORMITY

### Medical

---

### Declaration of Conformity

**Manufacturer:**  
**Address:**

**Product:** Disposable Nitrile Examination Gloves  
**Designation:** X-Small, Small, Medium, Large, X-Large, XX-Large

We herewith declare that the above-mentioned devices comply with the European Medical Device Regulation (EU) MDR 2017/745 and PPE Regulation (EU) 2016/425. The EU declaration of conformity is issued under the sole responsibility of the manufacturer:

By formulating the products, the chemical substances selected was nitriles, and compliance to REACH, RoHS, Halogen-Free, SVHC (RI).

We strict following the standard of U.S. and EU, no DEHP, BBP, DBP, and DEHP is using in any vinyl products.

**STANDARDS**  
Standards Harmonized Standards applicable to this product are:  
EN455-1, EN455-2, EN455-3, EN455-4,  
EN374-1, EN374-2, EN374-3, EN374-5

**Signature:** 

**Date:** March 02, 2021

### Medical

### EU Declaration of Conformity

**Manufacturer:**  
**Address:**  
**SRN:**  
**European Representative:** Lotus NL B.V.  
koningin Julianaplein10, Te Verd, 2585AA The Hague, Netherlands  
**SRN:**  
**Product:** Disposable Medical nitrile exam glove  
X-Small, Small, Medium, Large, X-Large

**GMDN Code:**  
**UMDN Code:**  
**Basic UDI:**

**Classification (MDR, Annex VIII):** Class I, Rule 1.  
**Conformity Assessment Route:** EU DECLARATION OF CONFORMITY following the Annex II + Annex III + Article 19 of MDR (EU) 2017/745.

We herewith declare that the above mentioned de products meet the transposition into national law, the provisions of the following EU Regulation and Standards. All supporting documentations are retained under the premises of the manufacturer:  
is exclusively responsible for the declaration of conformity.

**General applicable regulations, directives:**  
Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC.  
Applied standards, common specification, guidance:  
EN 455-1:2020, EN 455-2:2015, EN 455-3:2015, EN 455-4:2009

**Signature:** \_\_\_\_\_  
**Date:** July 05, 2021

### Chemical

### EU Declaration of Conformity - PPE

**The manufacturer**

**Representative (EU, Switzerland):** CURADEN AG, Amlehnstrasse 22, 6010 Kriens, Switzerland

**Declares under his sole responsibility, that the PPE reference** : described hereafter:

**Protective nitrile Glove**  
**PPE to be used against Category III risks**

EN ISO 374-1:2016-Type C

EN ISO 374-5:2016

is in conformity with the provisions of Regulation (EU)2016425 and with the European harmonized standards EN ISO 374-1+A1:2018, EN ISO 374-5:2016\*, EN ISO 21420:2020\* (\*see SATRA test report) and is identical to the PPE which is subject to the EU-Type examination, under certificate number 0197, issued by the Notified Body:  
TUV Rheinland LGA Products GmbH  
Tillystrabe 2  
90431 Nürnberg  
Germany

and is subject to the procedure set out in Annex VII (Module C2) of the Regulation under the supervision of the Notified body:  
TUV Rheinland LGA Products GmbH  
Tillystrabe 2  
90431 Nürnberg  
Germany

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Meditec Germany GmbH